

SELECTIVE “ILY” Type

Bone marrow aspiration needle

*Whistle tip needle, anatomic handle.
Removable stylet with whistle point perfectly
coupled at the cannula tip.*



*Punta a becco di flauto
Whistle tip needle*

	Ø Gauge	Ø mm	Lunghezza Length 50 mm	Lunghezza Length 70 mm	Lunghezza Length 100 mm	Lunghezza Length 120 mm	Lunghezza Length 150 mm
○	11	3,20	IT11-05	IT11-07	IT11-10	IT11-12	IT11-15
●	13	2,50	IT13-05	IT13-07	IT13-10	IT13-12	IT13-15
●	14	2,10	IT14-05	IT14-07	IT14-10	IT14-12	IT14-15
●	15	1,80	IT15-05	IT15-07	IT15-10	IT15-12	IT15-15
○	16	1,60	IT16-05	IT16-07	IT16-10	IT16-12	IT16-15
●	18	1,20	IT18-05	IT18-07	IT18-10	IT18-12	IT18-15

Box 10 units / Different sizes are available on request

SELECTIVE E “ILY” Type

Bone marrow explantation needle

Whistle tip needle with large side holes nearby the needle point, anatomic handle with Luer Lock connection.

Removable stylet with whistle point perfectly coupled at the cannula tip.



2 HOLES



Whistle tip needle

2 holes cannula

	Ø Gauge	Ø mm	Lunghezza Length 70 mm	Lunghezza Length 100 mm	Lunghezza Length 150 mm
○	11	3,20	ITX11-07	ITX11-10	ITX11-15
●	13	2,50	ITX13-07	ITX13-10	ITX13-15
●	14	2,10	ITX14-07	ITX14-10	ITX14-15
●	15	1,80	ITX15-07	ITX15-10	ITX15-15
○	16	1,60	ITX16-07	ITX16-10	ITX16-15
●	18	1,20	ITX18-07	ITX18-10	ITX18-15

Box 10 units / Different sizes are available on request