

EVO CUT

Automatic tru-cut biopsy system

Cutting needle tip. Automatic Tru-cut type system with, 20 mm depth penetrations. Depth marks on the cannula. Stylet connected to the system with notch to keep the specimen in, whistle and penetrating point.

*Device complete of **simple shot** with cannula retraction or **complete shot** with cannula and stylet retraction.*

Echogenic marker: inner



Whistle tip needle

10mm/20 mm notch

	Ø Gauge	Ø mm	Lunghezza Length 100 mm	Lunghezza Length 150 mm	Lunghezza Length 200 mm	Lunghezza Length 250 mm
●	14	2,10	EC14-10	EC14-15	EC14-20	EC14-25
○	16	1,60	EC16-10	EC16-15	EC16-20	EC16-25
●	18	1,27	EC18-10	EC18-15	EC18-20	EC18-25
●	20	0,90	EC20-10	EC20-15	EC20-20	EC20-25

Box 10 units / Different sizes are available on request