

IVS/IVH

GENERAL PURPOSE ONE/SHOT TYPE TEFLON CATHETERS

Needle with cannula for abscess drainage, cysts emptying, metal guide-wires positioning, special procedures.

IVH: needle with cannula (three pieces), Teflon catheter with 2 holes (final and side-holes) complete with needle and stylet

IVS: needle with cannula (three pieces), clear Teflon catheter with final hole complete with needle and stylet

DIRECT TECHNIQUE

IVH

CODE	Diameter	Length	Box
IVH 0515	Fr. 5.0(19G)	150mm	10
IVH 0520	Fr. 5.0(19G)	200mm	10
IVH 0525	Fr. 5.0(19G)	250mm	10
IVH 0615	Fr. 6.0(18G)	150mm	10
IVH 0620	Fr. 6.0(18G)	200mm	10
IVH 0625	Fr. 6.0(18G)	250mm	10
IVH 0630	Fr. 6.0(18G)	300mm	10
IVH 0715	Fr. 7.0(16G)	150mm	10

IVS

CODE	Diameter	Length	Box
IVS 0415	Fr. 4.0(23G)	150mm	10
IVS 0420	Fr. 4.0(23G)	200mm	10
IVS 0515	Fr. 5.0(19G)	150mm	10
IVS 0520	Fr. 5.0(19G)	200mm	10
IVS 0525	Fr. 5.0(19G)	250mm	10
IVS 0530	Fr. 5.0(19G)	300mm	10
IVS 0615	Fr. 6.0(18G)	150mm	10
IVS 0620	Fr. 6.0(18G)	200mm	10
IVS 0625	Fr. 6.0(18G)	250mm	10
IVS 0715	Fr. 7.0(16G)	150mm	10
IVS 0720	Fr. 7.0(16G)	200mm	10
IVS 0725	Fr. 7.0(16G)	250mm	10
IVS 0815	Fr. 8.0(14G)	150mm	10
IVS 0820	Fr. 8.0(14G)	200mm	10



PERCUTANEOUS DRAINAGES

M.S. Hospital Service S.p.A. is certified ISO 9001 and EN 46001