

Agile performance



Enhanced crossability

Optimum track and push

Fast deflation time

NIMBLE 

PTA OTW BALLOON - .035" WIRE COMPATIBLE

General features

- Over the wire balloon catheter (OTW)
- Usable Catheter Length (UCL): 80, 130 and 150 cm
- Sheath compatibility: from 5 to 7F
- Guidewire compatibility: .035"



Shaft technical data

- Material: polyamide
- Proximal shaft: dual lumen with large crescent shape
- Highly lubricous hydrophilic coating (SlipSkin™)

Balloon technical data

- Soft distal tip
- Material: Polyn P™ Polyamide
- Two platinum markers
- Nominal Pressure (NP): 6 and 8 atm
- Rate Burst Pressure (RBP): 10, 12, 13 and 14 atm

Ordering information

Catheter Length 80 cm				
OTW .035"	Length (mm)			
Diameter (mm)	20	40	60	80
5	ICNIM050020S	ICNIM050040S	ICNIM050060S	ICNIM050080S
6	ICNIM060020S	ICNIM060040S	ICNIM060060S	ICNIM060080S
7	ICNIM070020S	ICNIM070040S	ICNIM070060S	ICNIM070080S
8	ICNIM080020S	ICNIM080040S	ICNIM080060S	ICNIM080080S
9	ICNIM090020S	ICNIM090040S	ICNIM090060S	ICNIM090080S
10	ICNIM100020S	ICNIM100040S	ICNIM100060S	ICNIM100080S
12	ICNIM120020S	ICNIM120040S	ICNIM120060S	ICNIM120080S

Catheter Length 130 cm				
OTW .035"	Length (mm)			
Diameter (mm)	20	40	60	80
5	ICNIM050020M	ICNIM050040M	ICNIM050060M	ICNIM050080M
6	ICNIM060020M	ICNIM060040M	ICNIM060060M	ICNIM060080M
7	ICNIM070020M	ICNIM070040M	ICNIM070060M	ICNIM070080M
8	ICNIM080020M	ICNIM080040M	ICNIM080060M	ICNIM080080M
9	ICNIM090020M	ICNIM090040M	ICNIM090060M	ICNIM090080M
10	ICNIM100020M	ICNIM100040M	ICNIM100060M	ICNIM100080M
12	ICNIM120020M	ICNIM120040M	ICNIM120060M	ICNIM120080M

Catheter Length 150 cm				
OTW .035"	Length (mm)			
Diameter (mm)	20	40	60	80
5	ICNIM050020L	ICNIM050040L	ICNIM050060L	ICNIM050080L
6	ICNIM060020L	ICNIM060040L	ICNIM060060L	ICNIM060080L
7	ICNIM070020L	ICNIM070040L	ICNIM070060L	ICNIM070080L
8	ICNIM080020L	ICNIM080040L	ICNIM080060L	ICNIM080080L
9	ICNIM090020L	ICNIM090040L	ICNIM090060L	ICNIM090080L
10	ICNIM100020L	ICNIM100040L	ICNIM100060L	ICNIM100080L
12	ICNIM120020L	ICNIM120040L	ICNIM120060L	ICNIM120080L

Please contact your local sales representative for product availability



INOMED HEALTH
INNOVATIVE MEDICAL DEVICES

PERIPHERAL VASCULAR

Agile performance



Enhanced crossability

Optimum track and push


Balanced performance

NIMBLE 
EXTENDED

PTA OTW BALLOON - .035" WIRE COMPATIBLE

General features

- Over the wire balloon catheter (OTW)
- Usable Catheter Length (UCL): 80, 130 and 150 cm
- Sheath compatibility: from 5 to 7F
- Guidewire compatibility: .035"



Shaft technical data

- Material: polyamide
- Proximal shaft: dual lumen with large crescent shape
- Highly lubricous hydrophilic coating (SlipSkin™)

Balloon technical data

- Soft distal tip
- Material: Polyn P™ Polyamide
- Two platinum markers
- Nominal Pressure (NP): 6 and 8 atm
- Rate Burst Pressure (RBP): 10, 12, 13 and 14 atm

Ordering information

Catheter Length 80 cm				
OTW .035"	Length (mm)			
Diameter (mm)	100	120	150	200
4	ICNIM040100S	ICNIM040120S	ICNIM040150S	ICNIM040200S
5	ICNIM050100S	ICNIM050120S	ICNIM050150S	ICNIM050200S
6	ICNIM060100S	ICNIM060120S	ICNIM060150S	ICNIM060200S
7	ICNIM070100S	ICNIM070120S	ICNIM070150S	ICNIM070200S
8	ICNIM080100S	ICNIM080120S	ICNIM080150S	ICNIM080200S

Catheter Length 130 cm				
OTW .035"	Length (mm)			
Diameter (mm)	100	120	150	200
4	ICNIM040100M	ICNIM040120M	ICNIM040150M	ICNIM040200M
5	ICNIM050100M	ICNIM050120M	ICNIM050150M	ICNIM050200M
6	ICNIM060100M	ICNIM060120M	ICNIM060150M	ICNIM060200M
7	ICNIM070100M	ICNIM070120M	ICNIM070150M	ICNIM070200M
8	ICNIM080100M	ICNIM080120M	ICNIM080150M	ICNIM080200M

Catheter Length 150 cm				
OTW .035"	Length (mm)			
Diameter (mm)	100	120	150	200
4	ICNIM040100L	ICNIM040120L	ICNIM040150L	ICNIM040200L
5	ICNIM050100L	ICNIM050120L	ICNIM050150L	ICNIM050200L
6	ICNIM060100L	ICNIM060120L	ICNIM060150L	ICNIM060200L
7	ICNIM070100L	ICNIM070120L	ICNIM070150L	ICNIM070200L
8	ICNIM080100L	ICNIM080120L	ICNIM080150L	ICNIM080200L

Please contact your local sales representative for product availability



2012/09 IC0701555626/A