

CATCHSYSTEM

Bone Marrow biopsy set with CatchSystem device

Bone Marrow Biopsy Needle with crown needle tip and anatomic handle.

It includes the CatchSystem cradle which allows to avoid the deflection technique procedure. It includes the extraction wire and safety device to get the biopsy core out.

Extractable trocar tip stylet.



Trocar tip

Crown tip cannula

CatchSystem cradle

	Ø Gauge	Ø mm	Lunghezza Length 100 mm	Lunghezza Length 150 mm
●	08	4,00	CS08-10	CS08-15
●	09	3,50	CS09-10	CS09-15
○	11	3,20	CS11-10	CS11-15
●	13	2,50	CS13-10	CS13-15

Box 10 units / Different sizes are available on request